



ACM20

| Automatic cryostat

Details Info

Product specifications :

- MPU (Microprocessor) Control the temperature and LCD display automatically, and control the function of the section and trimming semi-automatically.
- Two Compressor for freezing (import from Denmark).Five freezing place : two freezing tray , freezing specimen holder , freezing box , freezing knife.
- Beautiful Streamline design
- Specimen retraction function ,protects specimen from blade injury.
- Adopt advance import actuation system for advance of the specimen holder .
- Defrost Function
- Automatic retraction system for section .
- Automatic retraction system for specimen holder .
- Anti-roll system which made the slice flat .
- slice the specimen flat ,thin ,even and continuously.

Technical Date: __

- Temperature of Freezing Box: -5°C~-35°C (±2°C)
- Temperature of Freezing Tray: -10°C~-40°C
- Temperature of Freezing specimen holder : -5°C~-40°C (±2°C)
- Temperature of Freezing Knife: -5°C~-25°C
- Range of Section Thickness
- Trimming thickness: 1-60um
- Setting value : from 1-60um 1um increment
- Maximum Sample: 55mm
- Maximum vertical specimen stroke:59mm
- Maximum horizontal specimen stroke :28mm via step motor
- Range of Angle Adjustment of freezing Knife: 0-15°
- Dimension: 600×640×1150mm
- Weight: 120kg
- Power : 220V/110V



Head Office :1600 AVENUE DE LORIMIER BUREAU, 384 MONTRÉAL (QUÉBEC) CANADA H2K3W8 USAOffice:128 W SUFFLOK AVE,CENTRAL ISLIP,11722 N
Middle east & Africa Office : warehouse 9 /shed 9, Technology park P.O. Box 54806 RAK FZE , UAE



Email Us At:

sales@acculabusa.com middleeast@acculabusa.com



Phone Support:

© Copyright 2014 by ACCULAB_USA - All rights reserved