

ACR-100

| Compressor

☰ Details Info

**Positive & negative pressure for medical use
Practicality and Merits of Oil-free Compressor for Medical Use**

Features:

- 1. Be able to supply oil-free, odorless, clean, and dry green gas;

Owing to its absence of the usage of oil, the product can generate healthy gas, which is beneficial to both patients and environment

- 2. Be able to be maintained with considerable convenience;

Without the using of oil, the product has been privileged in its easy and convenient daily maintenance, comparing with the other compressors

- 3. Be privileged in oral medical treatment;

The existence of oil molecules in the gas generated by the traditional compressor is likely to give rise to the deduction of medical effect. Conversely, the oil-free compressor is able to provide oil-free gas.


- 4. Be privileged in prolonging the service life;

Likewise, the oil molecules existed in the gas is likely to increase the emission of large amount of carbon deposition, which will in return accelerate the fraction of .the hand piece. However, the design of the oil-free compressor could avoid the above-mentioned deficiencybearing in

Specification:


Model	ACR-100
Voltage Volt/Hz	220/50/1

Head Model	Home-made high quality Head
Volume Flow	100L/min
0.4Mpa Production Rate	60L/min
Maximum Pressure	0.8Mpa
Pressure Switch Adjustment	0.55-0.8Mpa
Power	0.55KW
Rotational Speed	1400RPM
Container Capacity	24L
Net	24kg
Gross	32kg
Packing size (LWG)	500x500 x600m
Packing materiel	Cardboard Box
Fitting Dental Chair	1units

 Head Office :1600 AVENUE DE LORIMIER BUREAU, 384 MONTRÉAL (QUÉBEC) CANADA H2K3W8 USAOffice:128 W SUFFLOK AVE,CENTRAL ISLIP,11722 N Middle east & Africa Office : warehouse 9 /shed 9, Technology park P.O. Box 54806 RAK FZE , UAE

 Email Us At:

sales@acculabusa.com middleeast@acculabusa.com

 Phone Support:

© Copyright 2014 by ACCULAB_USA - All rights reserved