



ACR-400 / ACR-500 / ACR-600

| Compressor

☰ Details Info

Positive & negative pressure for medical use Practicality and Merits of Oil-free Compressor for Medical Use

Features::

1. Be able to supply oil-free, odorless, clean, and dry green gas;

Owing to its absence of the usage of oil, the product can generate healthy gas, which is beneficial to both patients and environment

2. Be able to be maintained with considerable convenience;

Without the using of oil, the product has been privileged in its easy and convenient daily maintenance, comparing with the other compressors

3. Be privileged in oral medical treatment;


The existence of oil molecules in the gas generated by the traditional compressor is likely to give rise to the deduction of medical effect. Conversely, the oil-free compressor is able to provide oil-free gas.

4. Be privileged in prolonging the service life;

Likewise, the oil molecules existed in the gas is likely to increase the emission of large amount of carbon deposition, which will in return accelerate the fraction of bearing in the hand piece. However, the design of the oil-free compressor could avoid the above-mentioned deficiency.


Specification:

Model	ACR-400	ACR-500	ACR-600
Voltage Volt/Hz	220/50/1		
Head	Fiat Head		
Volume Flow	300L/min		
0.4Mpa Production Rate	224L/min		
Maximum Pressure	0.8Mpa		
Pressure Switch Adjustment	0.55-0.8Mpa		
Power	2x1.1KW		
Rotational Speed	1400RPM		
Container Capacity	100L		
Net	95kg		
Gross	130kg		
Packing size (LWG)	1200x550 x1000mm		
Packing materiel	Wooden Box		
Fitting Dental Chair	4Units	5units	6units


 Head Office :1600 AVENUE DE LORIMIER BUREAU, 384 MONTRÉAL (QUÉBEC) CANADA H2K3W8 USAOffice:128 W SUFFLOK AVE,CENTRAL ISLIP,11722 N Middle east & Africa Office : warehouse 9 /shed 9, Technology park P.O. Box 54806 RAK FZE , UAE


 Email Us At:

sales@acculabusa.com
middleeast@acculabusa.com


 Phone Support:

© Copyright 2014 by ACCULAB_USA - All rights reserved