



ACR-700 / ACR-900

| Compressor

☰ Details Info

Positive & negative pressure for medical use Practicality and Merits of Oil-free Compressor for Medical Use

Features::

1. Be able to supply oil-free, odorless, clean, and dry green gas;

Owing to its absence of the usage of oil, the product can generate healthy gas, which is beneficial to both patients and environment

2. Be able to be maintained with considerable convenience;

Without the using of oil, the product has been privileged in its easy and convenient daily maintenance, comparing with the other compressors

3. Be privileged in oral medical treatment;

The existence of oil molecules in the gas generated by the traditional compressor is likely to give rise to the deduction of medical effect. Conversely, the oil-free compressor is able to provide oil-free gas.


4. Be privileged in prolonging the service life;

Likewise, the oil molecules existed in the gas is likely to increase the emission of large amount of carbon deposition, which will in return accelerate the fraction of bearing in the hand piece. However, the design of the oil-free compressor could avoid the above-mentioned deficiency.

Specification:

Model	ACR-700	ACR-900
-------	---------	---------

Voltage Volt/Hz	220/50/1	
Head	Fiat Head	
Volume Flow	400L/min	
0.4Mpa Production Rate	270L/min	
Maximum Pressure	0.8Mpa	
Pressure Switch Adjustment	0.55-0.8Mpa	
Power	2x1.5KW	
Rotational Speed	1450RPM	
Container Capacity	100L	
Net	111kg	
Gross	146kg	
Packing size (LWG)	1200x580 x1000mm	
Packing materiel	Wooden Box	
Fitting Dental Chair	7units	9units


 Head Office :1600 AVENUE DE LORIMIER BUREAU, 384 MONTRÉAL (QUÉBEC) CANADA H2K3W8 USAOffice:128 W SUFFLOK AVE,CENTRAL ISLIP,11722 N Middle east & Africa Office : warehouse 9 /shed 9, Technology park P.O. Box 54806 RAK FZE , UAE



Email Us At:

sales@acculabusa.com middleeast@acculabusa.com



Phone Support:

© Copyright 2014 by ACCULAB_USA - All rights reserved