



AHS62 - Automatic Immuno-Histo Stainer

| STAINER

☰ Details Info

Features:

Program can be changed at any time required and can be made according to the user's requirements.

More than one hundred programs can be made for standardized immunostaining and rapid special staining.

One computer can control several immuno-Histo Stainers simultaneously.

To ensure the continuity and completeness of an experimental run, an external power supply can be provided.

Such running conditions as running time and running progress can be clearly monitored at any time.

The standardization of the test method of immunostaining is becoming more and more important with the wide use of immunostaining technology in the modern pathologic diagnosis because it has a bearing on the prognosis and treatment of the patient.

AHS62 immunohistochemistry stainer adopts European advanced technology with an intelligent program interface which is easy to understand and operate, thus providing a scientific and reliable platform for the standardization of the test method of immunostaining.

Open the system, suitable for reagents of different

Suitable for various kinds of specimens; paraffin, refrigerated and punctured microtome sections, cell smear and marrow microtome section.

Precise micro dropping rates of reagent adjustable to the minimum of 50ul, thus saving expensive reagent.

Different dropping rates of antibody can be set for different tissue sections.

Excellent overall sealing keeps the inside environment wet.

Contents are under the protection of secret code.

The dropper point surface is Teflon treated and can be thoroughly cleaned to avoid cross contamination.

Movable scavenging gentle air blows off the liquid residue on the tissue sections.

Technical parameters:

Capacity	48 microtom section and 64 kinds of reagents
Reagent titer increment	50 μ l ,75 μ l ,100 μ l, 150 μ l, 200 μ l, 250 μ l
Outside dimensions(WXDXH)	1015X665X520 mm
Power supply	\pm 10% at 50/60HZ or 110V \pm 10% at 50/60HZ
Nominal power	350W
Weight	About 116 kg

**12x10ml****6x500ml****148x7ml****4x500ml****4x250ml****2x2x96holes****8x50ml****4x100ml****36x10ml**

Head Office :1600 AVENUE DE LORIMIER BUREAU, 384 MONTRÉAL (QUÉBEC) CANADA H2K3W8 USAOffice:128 W SUFFLOK AVE,CENTRAL ISLIP,11722 N
Middle east & Africa Office : warehouse 9 /shed 9, Technology park P.O. Box 54806 RAK FZE , UAE



Email Us At:

sales@acculabusa.com middleeast@acculabusa.com



Phone Support: