



## ALP-400

| Apex locator & pulp tester

☰ Details Info

### Feature:

1. Root canal apex finder and pulp vital measurement.
2. Wide color LCD screen, operation simple and friendly interface>
3. Fifth generation multi-frequency design for wet or dry accuracy, adjustment-free.
4. Progressive audio signal combined with on-screen graphics.
5. Rechargeable Li-ion Battery and charger.
6. Charge/working status intelligent switch to avoid EMI .

### Technical:

1. Adaptor voltage:AC100-240V,50/60HZ
2. Battery:7.4v/1000mAh changeable Li-ion battery.

3. Charge time:3hours

4. Display: color LCD screen ,buzzer.

**Volume/weight:**

volume(cm)	24x19x9
weight(KG/PCS)	1.5kg
Packing material	carton



Head Office :1600 AVENUE DE LORIMIER BUREAU, 384 MONTRÉAL (QUÉBEC) CANADA H2K3W8 USAOffice:128 W SUFFLOK AVE,CENTRAL ISLIP,11722 N Middle east & Africa Office : warehouse 9 /shed 9, Technology park P.O. Box 54806 RAK FZE , UAE



Email Us At:

sales@acculabusa.com middleeast@acculabusa.com



Phone Support:

---

© Copyright 2014 by ACCULAB\_USA - All rights reserved