



Battery Gel

| Horizontal Electrophoresis

☰ Details Info

Features

The ACBGEL battery mini-gel is designed specifically with safety and cost in mind: as no power supply is necessary. Results are obtainable within 1.5 hours if used with our completely safe educational electrophoresis kit, which is free of ethidium bromide. It includes non-toxic electrophoretic dyes, spotting pipettes, agarose, buffer and classroom protocols.

Technical Specifications

Gel dimensions (w x l)	6 x 7.5cm
Unit dimensions (w x l x h)	17 x 11 x 3.5cm
Max. sample (capacity)	32 samples
Buffer volume	50ml
Combs available	4, 8, 12 and 16
:No. of samples Thicknesses	1mm

Ordering Information

ACBGEL	Battery operated horizontal gel 6 x 7.5cm, 2 x 8 sample combs , casting gates.
ACBGEL-KIT	Discontinued
ACBGEL-CG	Casting gates Pk/2 for ACBGEL
ACBGEL-BAT	PP3 batteries for ACBGEL Pk/2
ACBG-1-4	Comb 4 sample, 1mm thick
ACBG-1-8	Comb 8 sample, 1mm thick
ACBG-1-12	Comb 12 sample, 1mm thick
ACBG-1-16	Comb 16 sample, 1mm thick



Head Office :1600 AVENUE DE LORIMIER BUREAU, 384 MONTRÉAL (QUÉBEC) CANADA H2K3W8 USAOffice:128 W SUFFLOK AVE,CENTRAL ISLIP,11722 N Middle east & Africa Office : warehouse 9 /shed 9, Technology park P.O. Box 54806 RAK FZE , UAE



Email Us At:

sales@acculabusa.com middleeast@acculabusa.com



Phone Support: