



CE4-CYT - Cell smear centrifuge

| Cell smear centrifuge

☰ Details Info

CE4-CYT cell smear centrifuge is the special equipment of liquid-based cell production.

Feature

1. All-steel construction with stainless steel centrifugal chamber.
2. To drive directly by large torque DC motor, speed range is 100 4000rpm, microcomputer control and simple operation.
3. The protection from overspeed, imbalance and door-cover, can insure the safe operation
4. Parameters can be modified at any time during the running with-out stopping the machine.
5. Accuracy of centrifugal force: quickly start, prompt speed reaching and automatic braking.

Technical Data

Model	CE4-CYT
Max Speed	4000r/min
Max RCF	2300×g

Max Capacity	12×0.5ml
Timer Range	0~99min
Interference Suppression	En61326
Max. Noise	≤50dB(A)
Power Supply	AC220v 50Hz
Rated Power	90W
Dimension	275×390×300mm(L×W×H)
Weight	20kg

Rotor

Rotor No.	Max Speed	Max RCF	Capacity
Angle Rotor No.1	3000r/min	730×g	6×0.5ml
Angle Rotor No.2	3000r/min	2300×g	12×0.5ml



Head Office :1600 AVENUE DE LORIMIER BUREAU, 384 MONTRÉAL (QUÉBEC) CANADA H2K3W8 USAOffice:128 W SUFFLOK AVE,CENTRAL ISLIP,11722 N
Middle east & Africa Office : warehouse 9 /shed 9, Technology park P.O. Box 54806 RAK FZE , UAE



Email Us At:

sales@acculabusa.com middleeast@acculabusa.com



Phone Support:

© Copyright 2014 by ACCULAB_USA - All rights reserved