



DCU-3050

| Dental Unit

☰ Details Info

Features:

Install the Plastic Antifouling Cover in instrument tray, main water and air switch
 Install the filter cleaning system in Air negative pressure suction and water negativepressure salivary ejector
 Continnence exteriorization pure water bottle.
 Air pressure locking system for easy manipulation.
 Adopt LINAK motor system from Denmark for stable movement.
 The outside armrest could open for patient.
 The headrest could adjust and lock

Standard accessories	
3 way syringe	1pc
Lamp with column	1set
X film viewer	1set
Clean water system	1set
Automatic warm water& rinsing system	1set
Motor driven dental chair	1set
Clinic stool	1set

Technical parameters	
Power Supply	50Hz/60Hz、 AC 220V
Air source	Air pressure 550-800kPa
Air source	Water pressure 200-400kPa

Packing size	
Size (mm)	1490×1020×1120
) kg (N/W)	210
) kg (G/W)	300



Head Office :1600 AVENUE DE LORIMIER BUREAU, 384 MONTRÉAL (QUÉBEC) CANADA H2K3W8 USAOffice:128 W SUFFLOK AVE,CENTRAL ISLIP,11722 N Middle east & Africa Office : warehouse 9 /shed 9, Technology park P.O. Box 54806 RAK FZE , UAE



Email Us At:

sales@acculabusa.com middleeast@acculabusa.com



Phone Support:

© Copyright 2014 by ACCULAB_USA - All rights reserved