



## DCU-3080

| Dental Unit

☰ Details Info

### Features:

- ◆ DCU-3080 the treating system composed of Bien-Air hand piece offers the best accuracy.
- ◆ Dental chair driver by import low noise DC motor results in smooth movement, and with spitting, reset and three preset positions.
- ◆ The lowest position is 410mm. Both the control box and the assistant console can be 90° rotated which creates more spare freedom.
- ◆ New integral ceramic spittoon, instrument tray and main chest are easy for cleaning; foot control has controlling system on instrument tray and assistant holder separately

Standard accessories	
3 way syringe	2pcs
Lamp with column	1set
X film viewer	1set
Clean water system	1set

Automatic warm water& rinsing system	1set
Motor driven dental chair	1set
Clinic stool	1set

Protecting device is installed on the backrest in case medical negligence

°rotated, and fixed time for flushing water Smooth glazed surface and ceramic spittoon provides more convenience for cleaning, and spittoon can be 90◆

◆Assistant console have a holder left for easy upgrade, dental chair ,rinsing &warm water can be controlled by touchpad.

◆ Instrument tray adopts air locking system which is easy for manipulation

Packing size	up	Down
Size (mm)	1680×1120×1120	
) kg (N/W	220	220
) kg (G/W	320	350

Head Office :1600 AVENUE DE LORIMIER BUREAU, 384 MONTRÉAL (QUÉBEC) CANADA H2K3W8 USAOffice:128 W SUFFLOK AVE,CENTRAL ISLIP,11722 N Middle east & Africa Office : warehouse 9 /shed 9, Technology park P.O. Box 54806 RAK FZE , UAE



Email Us At:

sales@acculabusa.com middleeast@acculabusa.com



Phone Support:

---

© Copyright 2014 by ACCULAB\_USA - All rights reserved