



DCU-3100

| Dental Unit

Details Info

Features:

The product is easy to be cleaned, which in turn could be guard against the potential interactional infection.

Apparatus board, assistant shelf, and oral lights could be shifted at both sides of dental chair with ease, thus bestly satisfying the dentist's needs in operation.

Standard accessories	
Fiber optic high speed handpiece Connector	1set
Fiber optic high speed handpiece Connector	1set
Fiber optic low speed handpiece Connector	1set

scalar	1set
curing light	1set
3 way syringe	2pcs
suction ejector	1set
saliva ejector	1set
clean water system	1set
LED operating light	1set
option X-film viewer	

Notes:

◆ and it is easy to be cleared. The suction apparatus adopts a concentrated pumping system, which warrants a high efficiency and convenient maintenance.

The assistant seat could be rotated in a large scope,

◆ function, reliable control, which could also be equipped at both the right and left side of the dental chair alternatively. Hand piece could be supplied water originated both from purified and pipe water resource. The operating board is characterized by its abundant

◆ The LED operating light is characterized by flexible roatation, downy light and sensor control, weak light source and steeples adjustment

◆ Spacious leather cushion, presenting its gracefulness and elegance. The newly designed power system could make the chair retain a minimum lower position possible. Moreover, it could also be rotated at a scope of 30 degree both clockwise and counter-clockwise. The dental chair is covered with comfortable and counter-clockwise.

◆Foot switch is able to control the motion of dental chair, oral lights, and the operating mode of hand piece.



Head Office :1600 AVENUE DE LORIMIER BUREAU, 384 MONTRÉAL (QUÉBEC) CANADA H2K3W8 USA Office:128 W SUFFLOK AVE,CENTRAL ISLIP,11722 N Middle east & Africa Office : warehouse 9 /shed 9, Technology park P.O. Box 54806 RAK FZE , UAE



Email Us At:

sales@acculabusa.com middleeast@acculabusa.com



Phone Support: