



DCU-3110

| Dental Unit

Details Info

Features:

DCU-3110 dental unit has up and down type. Main chest and instrument tray install Plastic Antifouling Cover and adopt the anticorrosion anti-rust, and fastness material.

- ◆ Install the filter cleaning system in Air negative pressure suction and water negative pressure salivary ejector.
- ◆ Install pure water bottle and lacking water alarm system.
- ◆ Adopt LINAK motor system from Denmark for stable movement. Dental unit size 410mm-750mm,weight 135kg.
- ◆ Install dental unit locking system.
- ◆ The dental backrest and the dental chair are equipped with Safety protection devices:when a chair and chair encountered obstacles in the process of decline,

the action is terminated and automatic. marginally to ensure safe, avoid accidents. S2311 dental unit has up and down type. Main chest and instrument tray install Plastic

Standard accessories	
Optical high speed handpiece	1pc
3 way syringe	2pcs
Lamp with column	1set

X film viewer	1set
Clean water system	1set
Automatic warm water& rinsing system	1set
Dental chair	1set
Doctor stool	1set

Packing size	Up	Down
Size (mm)	600×620×1090	1680×1120×1120
kg (N/W)	220	220
kg (G/W)	330	330

## Notes:

- ◆Dental unit is fitted with microcomputer controlling system, the main board could control programmable working positions, lamp, rinsing and flushing water, heating water
- ◆Operation light adopts cold light design, and less illumination surface, and more shadow less space can be acquired
- ◆ Only 9 seconds for backrest to move from lowest point to highest point
- ◆ The operation light would be controlled by Multi-functional foot control.
- ◆ 90° rotation of spittoon and assistant stand
- ◆ The dental cushion and backrest could disassemble for cleaning.



Head Office :1600 AVENUE DE LORIMIER BUREAU, 384 MONTRÉAL (QUÉBEC) CANADA H2K3W8 USAOffice:128 W SUFFLOK AVE,CENTRAL ISLIP,11722 N Middle east & Africa Office : warehouse 9 /shed 9, Technology park P.O. Box 54806 RAK FZE , UAE



Email Us At:

sales@acculabusa.com middleeast@acculabusa.com



Phone Support: