



DCU-3160

| Dental Unit

☰ Details Info

Features:

- ◆ DCU-3160 design using the functionality and the new art-design ideas realized the man-machine integration of perfect operation concept
- ◆ Dental unit using the classical built-in land design, is divided into upper, lower, the main chest and instrument tray adopt advanced injection molding process, anti-corrosion, antirust and fadeless.
- ◆ Air negative pressure aspirator and water pressure saliva ejector are equipped with cleaning mesh device, sediment filtering device.
- ◆ Dental chair power system adopts Denmark LINAK Linear mute motor, its work stability.
- ◆ The dental backrest and the dental chair are equipped with Safety protection devices: when a chair and chair encountered obstacles in the process of decline, the action is terminated and automatic marginally to ensure safe, avoid accidents
- ◆ Safety locking system between unit and chair.
- ◆ The dental chair is two armrests design, and the external armrest can be opened as requirements.
- ◆ Inner pure water system, and fulfill the requirements of the water-off and piping disinfection.
- ◆ The special large design of operation table could supply best working condition.
- ◆ Rotatable and anti-drop off hand piece holders fulfill the clinical needs and avoid hand piece drop off.
- ◆ DCU-3160 dental unit is fitted with microcomputer controlling system; the main board could control programmable working positions, lamp, rinsing and flushing water, heating water.
 - ◆ Rotatable hand piece holders apply more space for doctor. The doctor would operate on 9 clock and 12 clock position. S2316 dental unit is fitted with micro computer

Pull the hose holder slightly, hand piece would hold in any position with small resilience. ◆

- ◆ More convenience and humanized design of hand piece holder.
- ◆ Compact design of one layer type of instrument tray.
- ◆ More satisfied touch screen of control panel and LEDX-film viewer

Standard accessories	
3 way syringe	2pcs
Lamp with column	1set
X film viewer	1set
Clean water system	1set
Automatic warm water& rinsing system	1set
Motor driven dental chair	1set
Clinic stool	1set

Packing size	Up	Down
Size (mm)	1850×1180×1120	1700×1210×1120
) kg (N/W	225	225
) kg (G/W	360	360

Notes:

Adopt tow connectors, adjusted angle and reliable locking system which convenient for adults and children.

New integral ceramic spittoon Rotated main chest and spittoon

Multi-functional foot control can control operation light.

Assistant console have a holder left for easy upgrade, dental chair ,rinsing &warm water can be controlled by touchpad.



Head Office :1600 AVENUE DE LORIMIER BUREAU, 384 MONTRÉAL (QUÉBEC) CANADA H2K3W8 USAOffice:128 W SUFFLOK AVE,CENTRAL ISLIP,11722 N Middle east & Africa Office : warehouse 9 /shed 9, Technology park P.O. Box 54806 RAK FZE , UAE



Email Us At:

sales@acculabusa.com middleeast@acculabusa.com



Phone Support: