



## DCU-3180

| Dental Unit

### ☰ Details Info

#### Features:

- ◆ There are up-linked and down-linked design for the Dental unit, the main box and instrument plate adopt the advanced injection molding process, and anti-corrosion, anti-rust, fadeless. The instrument tray is equipped with transparent silicone antifouling disk protection mat, with 135 degrees high temperature sterilization disinfection, avoiding cross infection.
- ◆ New instrument arm with two sets ball beating ensures turn smooth, accurate position, hidden joint design, equipped with locking devices.
- ◆ The Automatic depressurization and pressurization purify water system, and satisfy the need water and pipeline disinfection, and is equipped with water alarm prompt device.
- ◆ New Embedded connection box design for Dental chair, the movements are safe and smooth, and chair in the range of 410mm to 750mm, load is more than 135kg.
- ◆ Graceful techniques, fashionable figuration and specification make the dental unit integrated into the global market.
- ◆ Only 9 seconds for backrest to move from lowest point to highest point.
- ◆ Install dental chair protection system. When dental chair gets obstacle during the movement, the dental unit stop and rise up a little.
- ◆ Install dental unit locking system.
- ◆ Double deck instrument tray provide more spare freedom.
- ◆ New instrument arm with two sets ball beating ensures turn smooth, accurate position, hidden joint design, equipped with locking devices.
- ◆ Integrated holder for curing light and scalar
- ◆ The operation light would be controlled by Multi-functional foot control.

#### Notes:

- ◆ Integrated health ceramic spittoon, no sharp edges smooth design, 90 ° rotation, easy four-handed operation.

- ◆ Imported electromagnetic valve control the mouthwash device, heating and setting the time.
- ◆ Upholstery made of soft material with new coupling system for easy replacement.
- ◆ Assistance convenient operation.

Standard accessories	
Fibre optic handpiece	1pc
3 way syringe	2pcs
Lamp with column	1set
X film viewer	1set
Clean water system	1set
Automatic warm water& rinsing system	1set
Motor driven dental chair	1set
Clinic stool	1set

Packing size	Up	Down
Size (mm)	1680×1120×1120	1680×1120×1120
Net Weight (N/W)	230	220
Gross Weight (G/W)	345	330


 Head Office :1600 AVENUE DE LORIMIER BUREAU, 384 MONTRÉAL (QUÉBEC) CANADA H2K3W8 USAOffice:128 W SUFFLOK AVE,CENTRAL ISLIP,11722 N Middle east & Africa Office : warehouse 9 /shed 9, Technology park P.O. Box 54806 RAK FZE , UAE



Email Us At:

sales@acculabusa.com middleeast@acculabusa.com



Phone Support:

---

© Copyright 2014 by ACCULAB\_USA - All rights reserved