



## EDT-500

| Endodontic treatment

☰ Details Info

Endodontic treatment (with apex locator)


### Feature:

1. Original Switzerland electro-motor
2. Folding color wide LCD screen
3. Drive-1:1 contra angle is available
4. Rotation speed and torque can be programmed to 9 settings
5. Auto torque reverse

- 6. Dual-frequency apex locator
- 7. 5 working modes
- 8. Auto apical reverse
- 9. Auto apical slow down
- 10. Auto start and stop
- 11. High capacity Li-ion chargeable battery
- 12. Auto power off and memory function

**Technical:**

- 1.adaptor working voltage:AC100-240V,50/60hz
- 2.9 settings for adjustable rotational speed ranging from 125 to 625 rpm
- 3.adjustable torque, range:0.3-5.2N.CM(3-52MN.M)



**Volume/weight:**

volume(cm)	31x22.5x9.5
weight(KG/PCS)	2.5KG
Packing material	carton

Head Office :1600 AVENUE DE LORIMIER BUREAU, 384 MONTRÉAL (QUÉBEC) CANADA H2K3W8 USAOffice:128 W SUFFLOK AVE,CENTRAL ISLIP,11722 N Middle east & Africa Office : warehouse 9 /shed 9, Technology park P.O. Box 54806 RAK FZE , UAE

Email Us At:  
 sales@acculabusa.com middleeast@acculabusa.com

Phone Support: