



## EMC34

---

| Embedding center

☰ Details Info

### Features:

- The newest high technology and Button Panel Control,new type of heating material.
- The machine is made up of embedding machine, freezing part .the design can be flexibly selected between embedding machine and freezing part.
- Automated processor control system make it possible for presetting the start time and temperature of any day.
- special nano technical make it heat up quickly and equably.
- The LCD display can show you the temperatures of the paraffin bath, thermostatic bath and paraffin nozzle. So the user can operator it easily.
- inductive head lamp and metal switch ,user can touch the switch by tweezers or knives to control the paraffin.
- Automatic memory the presetting temperatures and start to work based on the time.
- large capacity crock can contain enough paraffin for use at one time.
- multilevel security protection system.
- There are 3 heating hole for forceps in two side of the paraffin nozzle .

### Specification:

- capacity of the wax bath : 6000 ml.
  - Melting paraffin prefabricate: 0 ~99°C
  - Freezing range temperature: 0 ~ -20°C
  - Complete machine consuming more electricity: Less than 1KVA
- ~~ Control mode: Auto, foot
- Environment temperature: 0°C~ 40°C
  - Time control: Opening and closing time of the machine is set in 24 hours optionally.
  - nominal voltage: AC220V±10%.
  - rated frequency:50±1Hz
  - Demension: 1180×660×470mm



Head Office :1600 AVENUE DE LORIMIER BUREAU, 384 MONTRÉAL (QUÉBEC) CANADA H2K3W8 USA Office:128 W SUFFLOK AVE,CENTRAL ISLIP,11722 N  
Middle east & Africa Office : warehouse 9 /shed 9, Technology park P.O. Box 54806 RAK FZE , UAE



Email Us At:

[sales@acculabusa.com](mailto:sales@acculabusa.com) [middleeast@acculabusa.com](mailto:middleeast@acculabusa.com)



Phone Support:

---

© Copyright 2014 by ACCULAB\_USA - All rights reserved