



EMC34M

| Embdding center

Details Info

Features:

- The newest high technology and Touching Screen Control, new type of heating material.
- Small Freezing Spot new to the working stage ,more confidence and high-efficiency for user to quick freeze paraffin cassettes
- Automated processor control system make it possible for presetting the start time and temperature of any day.
- special Nano technical make it heat up quickly and equably.
- The Touch Screen can show you the temperatures of the paraffin bath, thermostatic bath and paraffin nozzle. So the user can operator it easily.
- inductive head lamp and metal switch ,user can touch the switch by tweezers or knives to control the paraffin.
- Automatic memory the presetting temperatures and start to work based on the time.
- Large capacity Paraffin Bath (6 L)can contain enough paraffin for use at one time.
- Multilevel security protection system.
- There are 3 heating hole for forceps in two side of the paraffin nozzle .

Specification:

capacity of the wax bath	6000 ml.
Melting paraffin prefabricate	0 ~99℃
Freezing Spot temperature	0 ~ -10℃
Complete machine consuming more electricity Less	than 1KVA

Control mode:	Auto, foot
Environment temperature	0°C~ 40°C
Time control:	Opening and closing time of the machine is set in 24 hours optionally.
Nominal voltage	AC220V/110V±10%.
Rated frequency	50/60±1Hz



Head Office :1600 AVENUE DE LORIMIER BUREAU, 384 MONTRÉAL (QUÉBEC) CANADA H2K3W8 USAOffice:128 W SUFFLOK AVE,CENTRAL ISLIP,11722 N
Middle east & Africa Office : warehouse 9 /shed 9, Technology park P.O. Box 54806 RAK FZE , UAE



Email Us At:

sales@acculabusa.com middleeast@acculabusa.com



Phone Support: