



EMC36

| Embedding center

☰ Details Info

Specification:

Model	EMC36
Nominal supply voltage	220V/50Hz
Maximum power draw	650VA
Working temperature	from ambient to 80°C increments 1°C
Capacity	Paraffin reservoir: 6.0 Lt Cassette/mold tray: 1.6 Lt
Dimensions	Width: 600mm Depth: 650mm Height: 460mm

	Weight: 30Kg
Cryo Console specification	Nominal supply voltage: 220V/50Hz Maximum power draw: 200VA Working temperature: from 0°C to -20°C increments 1°C
Dimensions	Width: 410mm Depth: 580mm Height: 440mm Weight: 25Kg
Thermal Console specification	Nominal supply voltage: 220V/50Hz Maximum power draw: 650VA Working temperature: from ambient to 80°C increments 1°C
Capacity	Paraffin tray: 5Lt
Dimensions	Width: 410mm Depth: 620mm Height: 440mm Weight: 13Kg



Head Office :1600 AVENUE DE LORIMIER BUREAU, 384 MONTRÉAL (QUÉBEC) CANADA H2K3W8 USA Office:128 W SUFFLOK AVE,CENTRAL ISLIP,11722 N
Middle east & Africa Office : warehouse 9 /shed 9, Technology park P.O. Box 54806 RAK FZE , UAE



Email Us At:

sales@acculabusa.com middleeast@acculabusa.com



Phone Support: