



EMC38C Cryo Console

| Embdding center

Details Info

Features

- The EMC 38-C Cryo Console consists an advanced compressor that can cool the paraffin blocks repidly . The Teflon coated surface guarantee a minimum of ice formation and easy clean. Temperature is adjustable from 15°Cto -15°C.

EMC 38-C Cryo Console can also be used as a stand alone unit near the microtome .

Specification :

Supply voltage	220V±10%, 50hz /110V±10%, 60hz
Power Draw	330W
Working temperature	from 15°C to -15 °C , increments 1°C
Capacity	70 cassettes
Dimension(W*D*H	390*540*390mm