



EMC38M Embedding Dispensing Console

| Embedding center

Details Info

Features:

1. Working surface :The tissue embedding cassettes and base mold trays are removable and can be set at different temperatures .
The ergonomic large working surface allows user to process more samples at the same time .The setting knobs offer more precise adjustment of the paraffin flow .
2. All functions on LCD control panle including :
Temperature of working surface
Temperature of molds, cassettes, warming trays
Temperature of paraffin reservoir
3. Weekly automatic On/Off programme .
4. The temperature can be set from ambient temperature to 90deg C and can be adjusted through the control panel
5. The large peltier cold spot is designed to fit both small and large samples . The temperature range is from room temp. to -5 degree .

Specification:

Nominal supply voltage	220V±10%, 50hz /110V±10%, 60hz
Nominal power	1200W
	Paraffin Tank: form ambient to 90°C

Working temperature (increments 1°C)	Cassettes/molds tray: from ambient to 90°C Working table : from ambient to 90°C Refrigerating table: from ambient to -5°C Forceps Holder: 65°C constant
Capacity	
Paraffin reservoir	6 Lt
Cassettes/molds tray	1.4 Lt(126 cassettes)
Dimensions Width Depth Height	560mm 640mm 400mm
Net Weight	About 26KGS



Head Office :1600 AVENUE DE LORIMIER BUREAU, 384 MONTRÉAL (QUÉBEC) CANADA H2K3W8 USA Office:128 W SUFFLOK AVE,CENTRAL ISLIP,11722 N
Middle east & Africa Office : warehouse 9 /shed 9, Technology park P.O. Box 54806 RAK FZE , UAE



Email Us At:

sales@acculabusa.com middleeast@acculabusa.com



Phone Support: