



FAM48 with separate control panel

| Microtome Fully-automatic

Details Info

Features:

- The casing of microtome is amde from fire-retsrdant ABS palstic with streamlined design .
- Spacious and ingenious waste tray.
- Advanced retraction function makes better sectioning. The retraction value is also adjustable .
- Universal and fast specimen holders fits international embedding cassettes.
- Handwheel locking may be activated at any position and improve operation safety standard .
- The electric coarse feed operates two speeds which is easy to record the feed.
- Separated control panel is much easier to operate.
- Three motorized sectioning modes (Cont, Single and Step) are available for user to select .
- Automatic memorized coarse feed and reload function can significantly improve cutting efficiency .

Specification:

	0 to 600μm
--	------------

Section thickness	0 to 2 $\mu$ m, 0.5 $\mu$ m increments 2 to 10 $\mu$ m, 1 $\mu$ m increment 10 to 20 $\mu$ m, 2 $\mu$ m increments 20 to 100 $\mu$ m, 5 $\mu$ m increments 100 to 600 $\mu$ m, 50 $\mu$ m increments
Trimming thickness	0 to 600 $\mu$ m 0.5 to 2 $\mu$ m, 0.5 $\mu$ m increments 2 to 10 $\mu$ m, 1 $\mu$ m increment 10 to 20 $\mu$ m, 2 $\mu$ m increments 20 to 100 $\mu$ m, 5 $\mu$ m increments 100 to 600 $\mu$ m, 50 $\mu$ m increments
Retraction thickness	5-100 $\mu$ m, 5 $\mu$ m increments (can deactivate)
Maximum specimen	40x50x30mm, or standard cassette
Specimen horizontal feed	28m
Specimen vertical feed	70mm
Specimen holder adjusted system	Horizontal orientation: $\pm 8^\circ$ Vertical orientation: $\pm 8^\circ$
Repositioning of knife holder base ,(left-right)	50mm
Electric Coarse Feed Speed	300 $\mu$ m/s (Slow) and 900 $\mu$ m/s (Quick)
Cutting speed	5 grades 1st: 100~164mm/s, 2nd: 164~228mm/s, 3rd: 228~292mm/s, 4th: 292~356mm/s, 5th: 356~420mm/s
Dimension:	Length: 550mm, Width: 410mm ,(without handwheel, 310mm) Height: 320mm ,
Weight	about 39 kg