



HGE12X12Q

| Horizontal Electrophoresis

☰ Details Info

Features:

The HGE5X8 is the smallest unit and is designed for low to medium numbers of samples. The small gel size is 48 x 75mm, and combs ranging from 3up to 9 samples. Simply by altering the combs, this compact unit is capable of resolving up to 9 different samples, prepping 1ml of sample or separating sample bands over a distance of 3cm. Injection molded construction - durable leak-proof environment for complete safety and long life.

Specification:

Gel Size	60×60mm; 120×60mm; 60×120mm 120×120mm..
Sample volume	3, 6, 8, 11, 13, 18 ,25.
Buffer Volume	550ml
Most can be performed simultaneously in 100 samples to electrophoresis	
The trays have the sample adding background color.	
Special construction effectively prevents tray drift	
All trays have handles.	

With the level adjustment.

Abundant buffer not only guarantees cooling effects, but also keeps the pH value stable during the whole experimental process. Unfolded electrodes make repair and parts replacement easy and quick.

The transparent top lid effectively prevents the liquid in the cell from volatilizing and avoids electricitm leakage. Auto-switch-off when the lid is removed

Head Office :1600 AVENUE DE LORIMIER BUREAU, 384 MONTRÉAL (QUÉBEC) CANADA H2K3W8 USAOffice:128 W SUFFLOK AVE,CENTRAL ISLIP,11722 N Middle east & Africa Office : warehouse 9 /shed 9, Technology park P.O. Box 54806 RAK FZE , UAE

Email Us At:

sales@acculabusa.com middleeast@acculabusa.com

Phone Support:

© Copyright 2014 by ACCULAB_USA - All rights reserved