



## HGE25X25

| Horizontal Electrophoresis

### Details Info

#### Features:

HGE25X25 is specially designed to meet most customer and market requirements. Because of sample quantity demands, our engineers completed several market surveys and receive feedback from many end users as part of our design process.

HGE25X25 is usually served for the currently electrophoresis available, features with the DNA electrophoresis screening. Several rows of samples for a total of 96 can be done simultaneously.

#### Specification:

Gel Size	245×250mm.
Sample volume	52
Buffer Volume	1700ml.
Abundant buffer not only guarantee cooling effects, but also keep the pH Value stable during the whole experimental process. Unfolded electrodes make repair and parts replacement easy and quick.	
The transparent top lid effectively prevents the liquid in the cell from volatilizing and avoids electricity leakage. Auto-switch-off when the lid is removed	