



## HGE5X8

| Horizontal Electrophoresis

☰ Details Info

### Features:

The HGE5X8 is the smallest unit and is designed for low to medium numbers of samples.

The small gel size is 48 x 75mm, and combs ranging from 3up to 9 samples. Simply by altering the combs,

this compact unit is capable of resolving up to 9 different samples,prepping 1ml of sample or separating sample bands over a distance of 3cm.Injection molded construction - durable leak-proof environment for complete safety and long life.

### Specification:

Gel Size	48×75mm.
Sample volume	3, 5, 9
Buffer Volume	150ml
Abundant buffer not only guarantees cooling effects, but also keeps the pH Value stable during the whole experiment process. Unfolded electrodes make repair and parts replacement easy and quick.	
The transparent top lid effectively prevents liquid in the cell from volatilizing and avoids electricity leakage.Auto-switch-off when the lid is removed.	