



MRM-31

| Microtome Rotary

☰ Details Info

Product specifications :

- New special and improved blade carrier which Non-stick with tissue section in surface
- High precision specimen feed
- Universal cassette clamp and standard specimen clamp
- The hand wheel can be locked in any position or highest position , The dual fail-safe control system is safe and reliable
- It adopt the special blade Carrier , it is more safe and convenient to use the disposable blade directly . good durability and long service life
- Smoothly running hand wheel
- Spacious, integrated section waste tray
- One New coming Special Blade carrier can use both of High and Lower Disposable Blade (Optional Accessories)

Technical Data:

Section thickness range	0- 60um
Setting values	From 0-2 um in 0.5 um –increments Form 2-10um in 1 um-increments Form 10-20 um in 2 um-increments Form 20-60um in 5 um-increments
Horizontal specimen stroke	mm 28

Vertical specimen stroke	60mm
Precision error	± 5%
Sample orientation	8° along the X-Y-axis, rotatable 360°
Maximum section size	50 × 45mm
Dimensions incl. tray (L/T/H)	300mm (L) X570mm (W) X270mm (H)
Net weight:	34kg



Head Office :1600 AVENUE DE LORIMIER BUREAU, 384 MONTRÉAL (QUÉBEC) CANADA H2K3W8 USA Office:128 W SUFFLOK AVE,CENTRAL ISLIP,11722 N
Middle east & Africa Office : warehouse 9 /shed 9, Technology park P.O. Box 54806 RAK FZE , UAE



Email Us At:

sales@acculabusa.com middleeast@acculabusa.com



Phone Support:

© Copyright 2014 by ACCULAB_USA - All rights reserved