



MRM-32

| Microtome Rotary

☰ Details Info

Features:

it don't need to lubricate and maintain the microtome frequently .

The hand wheel in the right can be locked in any position , and the handle can be locked in highest position , it is double safe for changing blade or tissue block .

The specimen Orienting system: axes X and Y 8° , axis 180°.

It adopt the special blade holder , it is more safe and convenient to use the disposable blade directly . it also can to longer the user life of disposable blade

Technical Data

| | |
|-----------------------|---|
| Slice thickness range | 0.5- 60um |
| Setting values | <p>from 0-2 um in 0.5 um – increments</p> <p>Form 2-10um in 1 um-increments</p> <p>Form 10-20 um in 2 um-increments</p> |

| | |
|----------------------------|---------------------------------|
| | Form 20-60um in 5 um-increments |
| Horizontal specimen stroke | 28 mm |
| Vertical specimen stroke | 60mm |
| Slip precision | ± 5% |
| Maximal slice section | 50 × 45mm |



Head Office :1600 AVENUE DE LORIMIER BUREAU, 384 MONTRÉAL (QUÉBEC) CANADA H2K3W8 USAOffice:128 W SUFFLOK AVE,CENTRAL ISLIP,11722 N
Middle east & Africa Office : warehouse 9 /shed 9, Technology park P.O. Box 54806 RAK FZE , UAE



Email Us At:

sales@acculabusa.com middleeast@acculabusa.com



Phone Support:

© Copyright 2014 by ACCULAB_USA - All rights reserved