



SCL-110

| scaler&curing light

Details Info

ultrasonic scaler

Features

The handpiece is made of the piezoelectric ceramic of ultrasonic.
 The handpiece has protective function when the output power is large
 The handpiece could be pulled and inserted freely. And it could be disinfected at 135°C,0.22Mpa high pressure.
 The adjustment and indication of the working strength in the same button. It is concise. The handpiece is compatible with Mectron

Specification

Power	100-240V~220V 50Hz or ~ 50/60Hz
Frequency	31kHz±3kHz
Weight	0.8kg
Dimension	180mm×150mm×90mm