



## SCL-135

| scaler&curing light

Details Info

ultrasonic scaler

### Features

The handpiece is made of the piezoelectric ceramic of ultrasonic.  
 The handpiece has protective function when the output power is large.  
 The handpiece could be pulled and inserted freely.And it could be disinfected at 135°C,0.22Mpa high pressure.  
 PC control, automatic frequency tracking, keeps the best work state Using 10 pcs LED to display the working strength  
 The hanpiece is compatible with EMS.

### Specification

Power	100-240V~220V 50Hz or ~ 50/60Hz
Frequency	29kHz±3kHz
Weight	1kg
Dimension	180mm×150mm×90mm