



SCL-80

| scaler&curing light

Details Info

ultrasonic scaler (built-in)

Features

Buit-in scaler,it is used for the dental chair.

The handpiece is made of the piezoelectric ceramic of ultrasonic.

The handpiece has protective function when the output power is large

The handpiece could be pulled and inserted freely. And it could be disinfected at 135°C,0.22Mpa high pressure.

PC control, automatic frequency tracking, keeps the best work state.

The handpiece is compatible with EMS.

Specification

Power	~24V 50/60Hz 1.3A
Operation frequency:	31kHz±3kHz
Weight	0.15kg
Dimension	63mm×50mm×33mm