



SWF20

| Tissue Bath & Dryer

☰ Details Info

Feature:

- Adopted high technology and Touching Screen Control, new type of heating material to ensure Small volume and light weight.
- It contains bath and baking processor and independent heat power without any affect each other .
- Automatical memorization and recover the last setting temperature stored to start when the power start again.
- If the power is switched off suddenly, the machine can go on work and correct mistake automatically when the power switch on again.
- high-tech control system make it precision in temperature and rapid heating and baking slice.

Technical specification:

Environment temperature	0~40°C
Power	≤550W
Power voltage	AC220V±10% 50HZ
Memory function	Automatically retain set temperature after operating
Temperature for roasting	indoor temperaturare-90°C(maximum) adjustable and automatically

size of bath plate	425×415×215mm 10kg
size of baking plate	425×415×215mm 10kg
Water Temperature for spreading: indoor temperature-90°C(maximum) adjustable and automatically thermostatic	



Head Office :1600 AVENUE DE LORIMIER BUREAU, 384 MONTRÉAL (QUÉBEC) CANADA H2K3W8 USAOffice:128 W SUFFLOK AVE,CENTRAL ISLIP,11722 N
Middle east & Africa Office : warehouse 9 /shed 9, Technology park P.O. Box 54806 RAK FZE , UAE



Email Us At:

sales@acculabusa.com middleeast@acculabusa.com



Phone Support:

© Copyright 2014 by ACCULAB_USA - All rights reserved