



TR - Series

| Ultraviolet products

☰ Details Info

Features:

Acculab transilluminators are suited for the visualization of fluorescent signals as well as for their documentation. The high transparent UV protective screen. Can open and close to any angle ,safe and convenient to operation. UV filter glass, quartz ultraviolet lamp, long service life, light uniform
Accomodate single, dual wavelength or UV/ White
Reduced electrical consumption
Unlimited filter life expectancy for 254nm, 312 nm and white light
Ondulex® reflector
High UV output
Stainless steel frame

Specifications :

Model	TR201	TR261-312	TR261-254	TR261-365	TR262	TR262M	TR262S	TR15UV/ WHITE
-------	-------	-----------	-----------	-----------	-------	--------	--------	------------------

Filter size mm	150x200	210x260	210x260	210x260	210X210	210X210	210X210	150x200UV/ 150X200 White light
Wave length nm	312	312	254	365	365 / 254	365 / 312	254 / 312	312/white light
No. of tube watt	4x8	6x15	6x15	6x15	6x15	6x15	6x15	6x15
Dual intenisty	-	High 100% & low 70%	High 100% & low 70%	High 100% & low 70%	High 100% & low 70%	High 100% & low 70%	High 100% & low 70%	High 100% & low 70%
Input Voltage	220-volt	220-volt	220-volt	220-volt	220- volt	220-volt	220-volt	220-volt
Frequency	50 Hz	50 Hz	50 Hz	50 Hz	50 Hz	50 Hz	50 Hz	50 Hz

Head Office :1600 AVENUE DE LORIMIER BUREAU, 384 MONTRÉAL (QUÉBEC) CANADA H2K3W8 USAOffice:128 W SUFFLOK AVE,CENTRAL ISLIP,11722 N Middle east & Africa Office : warehouse 9 /shed 9, Technology park P.O. Box 54806 RAK FZE , UAE

Email Us At:

sales@acculabusa.com middleeast@acculabusa.com

Phone Support:

© Copyright 2014 by ACCULAB_USA - All rights reserved