



USC120W – Series

| Ultrasonic Cleaner

☰ Details Info

Features:

- Digital display and setting ultrasonic cleaning time
- Display working countdown time
- Display actual temperature of the cleaning tank
- Temperature setting range: Ambient temperature-80°C
- Timing range: 1-99 minutes adjustable
- Equipped with a special stainless steel basket, reduce sound cover
- Shell and cover are made by stainless steel
- Equipment controlled by Single-chip controller
- Last working parameters can be stored in memory after power-off

Specification :

| | Inner tank | volume | Frequency | Power | Heating power | Adjustable temperature | Adjustabletime | Drainage | Rack |
|-------|------------|--------|-----------|-------|---------------|------------------------|----------------|----------|------|
| Model | | | | | | | | | |

| | L/W/H (mm) | L | KHz | W | W | °C | min | | √ |
|----------------|----------------|---|-----|-----|-----|------------|------|---|---|
| USC120W- 3L60C | 240*140*100 | 3 | 40 | 120 | 120 | Ambient-80 | 1-99 | / | √ |
| USC120W-5L60C | 240*140*150 | 5 | 40 | 120 | 120 | Ambient-80 | 1-99 | / | √ |
| USC180W-6L60C | 300*155*150 | 6 | 40 | 180 | 150 | Ambient-80 | 1-99 | / | √ |



Head Office :1600 AVENUE DE LORIMIER BUREAU, 384 MONTRÉAL (QUÉBEC) CANADA H2K3W8 USAOffice:128 W SUFFLOK AVE,CENTRAL ISLIP,11722 N Middle east & Africa Office
: warehouse 9 /shed 9, Technology park P.O. Box 54806 RAK FZE , UAE



Email Us At:

sales@acculabusa.com middleeast@acculabusa.com



Phone Support: