



## VGE-B Series

| Vertical- Electrophoresis

Details Info

### Features:

The installation of the glass plate in just 15 second , and does not need the knob, fast and convenient .

With the gel making devices modular , you can make multi slab gel at the same time

Safety cover open button design ,convenient opening the cover

Abundant buffer not only guarantees cooling effects, but also keeps the PH Value stable during the whole experiment process

Auto – switch –off when the lid removed

### Specification:

Model	VGE19X17B	VGE19X22B	VGE22X11B	VGE22X22B
Gel area (WXL)mm	160X150	160X205	186X95	186X205
Gel Plate area (WXL)mm	190X165	190X220	216X110	216X205
Gel thickness (mm)	1.0			
Amount of Gel	1~2*(piece)			
Buffer capacity ml	~1500	~3000	~1800	~3500
No.of samples	21~70		25,40,52	

Volume (LXWXH)mm	254X140X200	254X140X250	280X140X135	280X140X250
Weight	4.0(kg)	5.2(kg)	3.2(kg)	5.6(kg)



Head Office :1600 AVENUE DE LORIMIER BUREAU, 384 MONTRÉAL (QUÉBEC) CANADA H2K3W8 USAOffice:128 W SUFFLOK AVE,CENTRAL ISLIP,11722 N Middle east & Africa Office : warehouse 9 /shed 9, Technology park P.O. Box 54806 RAK FZE , UAE



Email Us At:

sales@acculabusa.com middleeast@acculabusa.com



Phone Support:

---

© Copyright 2014 by ACCULAB\_USA - All rights reserved