



## VGE10X10Q Mini-Protean Electrophoresis Cell

| Vertical- Electrophoresis

☰ Details Info

### Features

VGE10X10Q Vertical Protein Tank provides the opportunity to run one, two or up to four gels simultaneously by using a glass plate sandwich.

Users can choose using a casting base or making gels on the plates according to personal habits.

Under the same experimental conditions, operators can use one to four plates for PAGE electrophoresis experiments simultaneously.

Safety design electrode prongs ensure the experiment will run under safe conditions.

### Specification:

Glass plate size	100×100mm.
Gel size	82×88mm
Sample volume	11, 15
Buffer Volume	750ml
Auto-switch-off when the lid is removed	