



VGE22X20 Rapid SSR

| Vertical- Electrophoresis

☰ Details Info

Feature:

Abundant buffer not only guarantees cooling effects, but also keeps the pH Value stable during the whole experimental process. With special gel caster, and realize the integrated casting gel

, Avoid the damage of gel out of second movement. Cam fastening methods,

this is the way to make the process of installing glass plate becomes very simple.

Gel area optimization, sample design is reasonable, the electrophoresis tank is widely used in the SSR amplification products of polyacrylamide gel electrophoresis (PAGE) detection.

Switching off when opening the lid, safe electrical cables and made the operation simpler

Specification:

Glass plate size	214×200 mm..
Gel size	186×95 mm.
Sample volume	25, 40, 52
Buffer Volume	1800ml