



## WTB21

| Tissue Bath & Dryer

☰ Details Info

### Specification:

- Adopted high technology and nice screen control, new type of heating material to ensure Small volume and light weight.
- Automatic memory of setting temperature and running .
- High precision in temperature controlling. heating fast ,durable and energy-serving .
- If the power is switched off suddenly, the machine can go on work and correct mistake automatically when the power switch on again.
- The LCD display the current temperature and presetting temperature. and make the communication between user and machine easily.

### Technical Date:

Temperature prefabricate	0~100°C set arbitrarily
Rating Power	300W
Thermostatic control with an accuracy of $\pm 0.5^{\circ}\text{C}$	
Power voltage	AC220V $\pm 10\%$ 50HZ
Dimension	425×415×215mm

Weight	10k
--------	-----



Head Office :1600 AVENUE DE LORIMIER BUREAU, 384 MONTRÉAL (QUÉBEC) CANADA H2K3W8 USAOffice:128 W SUFFLOK AVE,CENTRAL ISLIP,11722 N  
Middle east & Africa Office : warehouse 9 /shed 9, Technology park P.O. Box 54806 RAK FZE , UAE



Email Us At:

[sales@acculabusa.com](mailto:sales@acculabusa.com) [middleeast@acculabusa.com](mailto:middleeast@acculabusa.com)



Phone Support:

---

© Copyright 2014 by ACCULAB\_USA - All rights reserved