



AC-HD fully automatic microcomputer

| Horizontal Autoclave

Details Info

Characteristics:

- 1.For 4~6 minutes rapidly sterilizing.
- 2.Digital display of working status,touch type key.
- 3.With 3 fixed cycles of water adding,temperature rising,sterilizing,drying steam discharging controlled automatically.
- 4.The whole process of water adding,temperature rising,sterilizing,drying steam discharging controlled automatically.
- 4.Steam-water inner circulation system:no steam discharge,and the environment for sterilizing will be clean and dry.
- 5.Exhaust the cool air automatically.
- 6.Safe protection of water lacking.
- 7.Door safety lock system.
- 8.With four stainless steel sterilizing plates.
- 9.The chamber of the sterilizer is made of stainless steel.
- 10.Automatically shut off with beep reminding after sterilization.

11.With drying function.

Technical data:

MODEL	AC35-HD	AC50-HD
Sterilizing chamber volume	35L (φ300×500 mm)	50L (φ340×550 mm)
Working pressure	0.22Mpa	
Working temperature	134°C	
Adjustment of temperature	105-134°C	
Timer	0-99 min	
Heat average	≤ ± 1°C	
Power	2.7 KW / AC220V 50Hz	
Sterilizing plate mm	400×200×30 (4piece)	500×250×30 (4piece)
Dimension mm	730×550×510	770×600×540
Package dimension mm	810×680×650	880×730×700
G. W/N.W	90/73 kg	105/86 kg



Head Office :1600 AVENUE DE LORIMIER BUREAU, 384 MONTRÉAL (QUÉBEC) CANADA H2K3W8 USAOffice:128 W SUFFLOK AVE,CENTRAL ISLIP,11722 N Middle east & Africa Office : warehouse 9 /shed 9, Technology park P.O. Box 54806 RAK FZE , UAE



Email Us At:

sales@acculabusa.com middleeast@acculabusa.com



Phone Support: