



AC-HT Desktop Rapid Steam

| Horizontal Autoclave

☰ Details Info

Features

- 1, the European B standard, with three pulsating pre-vacuum and vacuum drying. Residual moisture $\leq 0.2\%$.
- 2, the LCD operation, touch buttons.
- 3, vacuum up -0.8bar , applicable to all types of packaging, non-packaged, solid, A type hollow, porous and internal piping equipment sterilization; ensure efficient removal of hollow instruments (such as adjusting the tooth Drill phone) inside the air, so high-temperature steam can fully reach every corner.
- 4, full computer control, modular, easy to operate interface.
- 5, with BOWIE & DICK tests for the detection of steam penetration.
- 6, with the degree of vacuum test at any time to ensure the detection of vacuum performance.
7. Fast independent steam generator.
- 8, built-in mini-printer, a detailed record of the sterilization data.
- 9, with a safety valve, manual safety locks, pressure locking system self-test failure, automatic over-temperature overpressure protection.
- 10, can be inactivated HIV (HIV), hepatitis B (HBV), mad cow virus and Bacillus.

Specification:

Model	AC35-HT	AC50-HT
Sterilizing Chamber volume	35L Ø300X500mm	50L Ø340X550mm

Working pressure	0.22Mpa	
Working temperature	105-134 ⁰ C	
Timer	0-99min	
Heat average	≤±1 ⁰ C	
Power	AC220V50HZ/27KW	
Dimension(mm)	730X550X510	770X600X540
Exterior package dimension	810X680X650	880X730X700
G.W/N.W	115/93kg	130/106kg



Head Office :1600 AVENUE DE LORIMIER BUREAU, 384 MONTRÉAL (QUÉBEC) CANADA H2K3W8 USAOffice:128 W SUFFLOK AVE,CENTRAL ISLIP,11722 N Middle east & Africa Office : warehouse 9 /shed 9, Technology park P.O. Box 54806 RAK FZE , UAE



Email Us At:

sales@acculabusa.com middleeast@acculabusa.com



Phone Support: