



ACD-BD Series Digital control with printer

| Horizontal Autoclave

☰ Details Info

Specification:

Model	ACD-12BD	ACD-18BD	ACD-22BD
Description	Small Autoclave	Small Autoclave	Small Autoclave
Class	B	B	B
Standard	BS EN 13060:2004 + A2:2010		
Plug	3 hole plug	3 hole plug	3 hole plug
Wattage	1.7 KW	1.7 KW	1.7 KW
Sterilization temperature	121°C, 134°C		

Product size	66 X 44.5 X 39CM	66 X 44.5 X 39CM	70 X 44.5 X 39CM
Diameter of chamber	20 CM	25 CM	25 CM
Depth of chamber	36 CM	35 CM	45 CM
Volumn of chamber	12 L	18 L	22 L
Material of chamber	Stainless steel 304		
Min size of table for installation	46 X 32 CM	46 X 32 CM	48 X 32 CM
Screen type	LED	LED	LED
Screen size	2.5 x 2cm, 3.8x 2cm	2.5 x 2cm, 3.8x 2cm	2.5 x 2cm, 3.8x 2cm
Language available	English	English	English
Lock system of door	Mechanic + electrical lockers(optional)		
Output	Printer (Optional)	Printer (Optional)	Printer (Optional)

Funcation :

Automatic Sterilization

5 sterilization program

B&D / Helix / Vacuum testing program

Dry time adjustable: 1-60 min

Buil-in printer (optional)

High quality bacteria filter

Big water tank allows 7-8 sterilization programs

Program	121°C Solid	121°C Universal	134°C Solid	134°C Universal
Vacuum Times	1	3	1	3
Temp (°C)	121	121	134	134
Pressure(KPA)	110	110	210	210
Time(Min.)	20	20	4	4
Program	Cotton	BD/Helix Test	Vacuum Test	
Vacuum Times	3	3	3	
Temp (°C)	134	134	/	
Pressure(KPA)	210	210	-80	
Time(Min.)	4	3.5	15	