



## ACZ-4246 Four-channel

| Coagulation Analyzer

☰ Details Info

### Features:

Testing channels: 4channel; 4 different parameters analysis in the same time.  
Sample pre-warming position:24 positions;  
reagent pre-warming position:6 positions  
Advanced scattered light principle and percentage analysis ensure accurate results.  
Low reagent consumption less than 20uL, open reagent.  
Special testing cup location system, adding sample bracket  
Advanced optical light matching system ensures precision result.  
Optional PT Derived method for testing fibrinogen, decrease the dosage of reagent  
Top quality thermostat, precision 37 , ±0.3  
Dispense precision CV<2%  
Memory for 10.000 testing result  
Internal thermal-sensitive printer  
240\* 128 large LCD display.  
Perform all routine tests such as PT. APTT. TT Fibrinogen and clotting factors.  
RS232 interface, PC connecting