



AHA-3231

| Hematology Analyzer

☰ Details Info

Specification :

Low sample consumption	venous 9.8 ul, prediluted 20 ul for twice testing one time
Large storage capacity	up to 10,000 samples +3 histograms
Adanced valve technology,long life	
Volume measurement by time ,not wrong warning	
3-part differentiation of WBC, 23 parameters, single channel counter	
Internal thermal-sensitive printer or external printer	
Automatically sample probe cleanning (inside and outside)	
RS232 interface, PC connecting	

Automatic diluting , Mixing , rinsing and clog clearing
Windows operation system graphical buttons mouse and keyboard operation
8.4" color TFT, Windows interface all testing parameter displayed simultaneously
Electrical resistance for counting and SFT method of Hemaglobin

📍
Head Office :1600 AVENUE DE LORIMIER BUREAU, 384 MONTRÉAL (QUÉBEC) CANADA H2K3W8 USAOffice:128 W SUFFLOK AVE,CENTRAL ISLIP,11722 N Middle east & Africa Office : warehouse 9 /shed 9, Technology park P.O. Box 54806 RAK FZE , UAE

✉
Email Us At:

sales@acculabusa.com middleeast@acculabusa.com

☎
Phone Support:

© Copyright 2014 by ACCULAB_USA - All rights reserved