



BBR-B Series Single Door

| Blood Bank

☰ Details Info

Blood Bank Refrigerator

Blood Bank Refrigerator is specialized refrigerating equipment for blood (whole blood) storage application, which has been used in medical industries as well as cold storage of pharmaceuticals, biological products, vaccine, and so forth. It is suitable for hospitals, blood banks and sanitation and ant epidemic stations.

ADVANTAGE:

1. Four sensors design: two for display and two for controlling.
2. Foam-glass door, easy to observe ,also keep a better cooling temperature.
3. Forced air refrigeration system,ensure inner temperature constant and uniform.
4. Adjustable shelves, improve space utilization; blood storage baskets can be classified.
5. Individual transparent inner door, effectively reduce cool air loss when opening the door.

TECHNICAL PARAMETERS:

Model	BBR-120B	BBR -160B	BBR -210B	BBR -300B	BBR -400B
Capacity	120L	160L	210L	300L	400L
Temp. Range	4 ±1°C				
Temp. Accuracy	0.1°C				

Control System	Microprocessor control				
Display	Large screen LED display				
Alarm	Audio and visual alarm for: High or low temperature, Sensor failure, System failure; Door ajar; Power failure (all models) Condenser abnormal (for model:BBR-160B/210B) Low battery (for model: BBR-120B/300B/400B)				
Refrigeration Type	Forced air refrigeration system				
Refrigerant	R134a, CFC-Free				
Compressor	International famous compressor				
Condenser & Evaporator	High efficiency air-cooled condenser, made of copper				
Defrost	Auto defrost				
Construction	Unibody Design, Mono-assembly foaming (Rigid polyurethane insulation material)				
Internal Material	Stainless steel	Stainless steel	Stainless steel	Cold-rolled steel with antibacterial powder coating	Stainless steel
External Material	Cold-rolled steel with anti-bacteria powder coating				
Door	Three-ply transparent glass door with lock, Heating film preventing frosting, Sealing joint strip with magnetic force (for model: BBR-160B/210B)				
	Inner and outer double door design, hollow transparent glass outer door, transparent plastic inner door (for model: BBR-120B/300B/400B)				
Shelves	Adjustable shelves				
Casters	4 casters				
Consumption	220W	180W	200W	403W	430W
Power Supply	AC220V±10%, 50/60Hz				
External Size (W*D*H),mm	598*618*1110	522*600*1300	522*600*1550	660*644*1960	750*761*1799
Standard Accessory	2 shelves, 4 baskets, fluorescent light	Temperature printer, 3 shelves, 9 baskets, fluorescent	Temperature printer, 4 shelves, 12 baskets, fluorescent	5 shelves, 15 baskets, fluorescent light	5 shelves, 20 baskets, fluorescent light