



BBR-L Series

| Blood Bank

☰ Details Info

Main Characteristics:

- 1 Digital temperature recorder.
- 2 Large screen LED display
- 3 Audio and visual alarm
- 4 Three-ply transparent glass door with lock, Heating film preventing frosting, Stealing joint strip with magnetic force
- 5 Unibody design with round angel structure, easy for cleaning
- 6 Drawers/adjustable shelves made of quality steel wire are convenient to store articles.

Specification:

Model	BBR-120L	BBR-250L	BBR-300L	BBR-450L
Capacity	120L	250L	300L	450L

Temp Range	4 ±1°C			
Temp Precision	0.1°C			
Control System	Microprocessor control			
Display	Large screen LED display			
Alarm	Audio and visual alarm for: Height or low temperature; System failure; Door ajar; Power failure			
Refrigeration Type	Forced air refrigeration system			
Refrigerant	R134a, CFC-Free	R404a, CFC-Free	R134a, CFC-Free	
Compressor	International famous compressor			
Condenser & Evaporator	Height efficiency air-cooled condenser, made of copper			
Defrost	Auto defrost			
Construction	Unibody Design, Mono-assembly foaming(Rigid polyurethane insulation material)			
Internal Material	Stainless steel	ABS plate	Cold-rolled steel with anti bacteria powder coating	
External Material	Cold-rolled steel with anti-bacteria powder coating			
Door	Three-ply transparent glass door with lock Heating film preventing frosting Stealing joint strip with magnetic force			
Drawers/Shelves	Shelves	SS drawers	Adjustable shelves made of PVC coated wire	
Casters	4 casters			
Consumption	160W	350W	500W	450W
Power Supply	110/220V±10%, 50/60Hz			
External Size(L*W*H)	520*600*1050mm	570*530*1600mm	615*605*1780mm	660*700*1980mm
Standard Accessory	Temperature printer, 2 shelves baskets, fluorescent light	Digital temperature recorder, 6 drawers fluorescenlight	Temperature printer, 4 shelves, fluoresce light	5 shelves, fluorescent light

Head Office :1600 AVENUE DE LORIMIER BUREAU, 384 MONTRÉAL (QUÉBEC) CANADA H2K3W8 USAOffice:128 W SUFFLOK AVE,CENTRAL ISLIP,11722 N Middle east & Africa Office : warehouse 9 /shed 9, Technology park P.O. Box 54806 RAK FZE , UAE

Email Us At:

sales@acculabusa.com middleeast@acculabusa.com