



## BFA-100 Portable Body Fat Analyzer

---

| Small Equipments

☰ Details Info

### Specification

Mini Portable Handheld body fat analyzer

\*Body fat measurement

\*LCD display, easy to read

\*9 Person profile memory

### Technical Data

Display parameters: body fat rate, BMI, basal metabolism and body shape

BMI: 2.5 to 90.0

Basal metabolism: 385 kcal - 5000 kcal

Size: thin, standards, the hidden obesity, muscle and fat

Power supply: DC 3 v (2 section 7 battery)

Battery life: measurement completed about 2000 times

Dimensions: length: 210 MMX wide: 42 MMX high: 145 mm

Net weight: 260 g

Accessories: 2 section 7 battery, manual, save the card, certificate of approval

Applicable people: people love, heat by the movement of people, and the health-conscious people to lose weight

Temperature: + 10 degree centigrade to + 40 degree centigrade