



## BLT-180

| Blood Lead Test System

Details Info

### Features:

Blood Lead Test System is a simple, rapid, convenient and inexpensive blood lead detection system. The electrochemical method, differential potentiometric stripping analysis (DPSA) is used for quantitative determination of ppb level lead in human body fluid, including whole blood. Blood lead testing system is fast, with results in just 180 seconds. The system has stand-alone mode, and it can also work with the computer and printer to analyze the data and print the report. The system can save last 400 inspection records at most.



Head Office : 1600 AVENUE DE LORIMIER BUREAU, 384 MONTRÉAL (QUÉBEC) CANADA H2K3W8 USA Office: 128 W SUFFLOK AVE, CENTRAL ISLIP, 11722 N Middle east & Africa Office : warehouse 9 / shed 9, Technology park P.O. Box 54806 RAK FZE, UAE



Email Us At:

[sales@acculabusa.com](mailto:sales@acculabusa.com) [middleeast@acculabusa.com](mailto:middleeast@acculabusa.com)



Phone Support:

© Copyright 2014 by ACCULAB\_USA - All rights reserved