



CAE7X9

Clinical Electrophoresis

Details Info

Features:

The CAE7X9 is updated from The CAE7X9A with a double-bridge style design. This increased the samples two times compared with the CAE7X9A .

The Platinum filament are removable, this design make the replacement easily. The newly updated electric cables are more durable and adaptive to many kinds of power supply.

CAE7X9 is equipped with two sets of adjustable beams for controlling run distance.

Specification:

Bridge Size	1: 66×250mm×2 (W×L).
Bridge Size	2: 82×250mm×2 (W×L)
Buffer Volume	1000ml
Size of Cellulose Acetate Membrane	70×90