



ESR-10

| ESR Analyzer

Details Info

Features :

- 1.Insert&detect the sample at any time.
- 2.Display and print samples' No.(serial number).
- 3.Stored not less than 150 samples of the test results.
- 4.Power-off protection function.
- 5.Built-in thermal printer preparation equipment with automatic , manual printing options
- 6.Film keyboard operation.
- 7.RS-232 serial interface.

Specification :

Number of test samples	10 items at the same time
Max test speed	10samples/hour
Suited ESR tube	1.6ml
Sample volume	1.6ml

Read resolution	±0.2nm
Accuracy range of sample	50-64nm
Power	12V
Display	128×240 LCD screen
Accuracy range of sample	50nm~64nm

Head Office :1600 AVENUE DE LORIMIER BUREAU, 384 MONTRÉAL (QUÉBEC) CANADA H2K3W8 USAOffice:128 W SUFFLOK AVE,CENTRAL ISLIP,11722 N Middle east & Africa Office : warehouse 9 /shed 9, Technology park P.O. Box 54806 RAK FZE , UAE

Email Us At:
sales@acculabusa.com middleeast@acculabusa.com

Phone Support: