



ESR-20

| ESR Analyzer

Details Info

Features:

- ESR (westergren and wintrobe Value) and HCT
- Automatically active and automatic measuring item detection
- Random loading with vacuette blood collection tube compatible
- ESR real time result and curve display
- Dual data port supporting HIS/LIS
- Immediate, patch and add more printing supported ,the ESR curve printing available

Specifications:

Measuring Principle	By infrared Kinetic photometry
Test Channels	20
Measuring Items	ESR (wintergreen and winrobe Value) HCT
ESR Measuring Range	(0-160)mm/h

ESR Measuring Accuracy	±1mm
HCT Measuring Range	0.2~1
HCT Measuring Accuracy	±0.03
CV	≤3%
Sample Volume	1.5ml
ESR Measuring Time	30 minute
HCT Measuring Time	<1 minute
Date Inquiry	N/A
Loading pattern	Random
Active Pattern	Automatically
Input	Keypad
Data Storage	More than 255 (HCT and ESR)
Screen	LCD
Display	ESR real time result and curve on LCD
Printer	Internal thermal printer
Printing Content	Result& curve
Data port	RS-232 port supporting HIS/LIS
Data uploading	Manual



Head Office :1600 AVENUE DE LORIMIER BUREAU, 384 MONTRÉAL (QUÉBEC) CANADA H2K3W8 USAOffice:128 W SUFFLOK AVE,CENTRAL ISLIP,11722 N Middle east & Africa Office : warehouse 9 /shed 9, Technology park P.O. Box 54806 RAK FZE , UAE



Email Us At:

sales@acculabusa.com middleeast@acculabusa.com



Phone Support: