



## Electronic sphygmomanometer

---

| Small Equipments

☰ Details Info

Memory function: 60 groups  
Pressurization / decompression: full-automatic (intelligent)  
Exhaust mode: automatically and quickly exhaust  
accuracy :blood pressure:within  $\pm 3$ mmHg; pulse number: within  $\pm 5\%$   
measuring range: 0 ~ 299mmHg; pulse : 40 ~ 200 times / minute  
the wrist circumference: 135mm ~ 195mm

📍  
Head Office :1600 AVENUE DE LORIMIER BUREAU, 384 MONTRÉAL (QUÉBEC) CANADA H2K3W8 USAOffice:128 W SUFFLOK AVE,CENTRAL ISLIP,11722 N Middle east & Africa Office  
: warehouse 9 /shed 9, Technology park P.O. Box 54806 RAK FZE , UAE

✉  
Email Us At: