



HBM-12 Hemoglobin Meter

| Small Equipments

☰ Details Info

Features:

Hemoglobin Meter

HB Meter

Principle: Optical reflectance

Sample: Capillary or venous whole blood (15~20 ul)

Specifications:

Parameter	Hemoglobin
Principle	Optical reflectance
Sample)	Capillary or venous whole blood (15~20 ul)

Strip	H12 Hemoglobin Test Strip
Speed	Result in 10~30 seconds
Range	4.0g/dL~24.0g/dL
Memory	250 test results
Calibration	Factory Calibrated according to ICSH
Automatization	Self-checking, auto-judging and displaying malfunction

Head Office :1600 AVENUE DE LORIMIER BUREAU, 384 MONTRÉAL (QUÉBEC) CANADA H2K3W8 USAOffice:128 W SUFFLOK AVE,CENTRAL ISLIP,11722 N Middle east & Africa Office : warehouse 9 /shed 9, Technology park P.O. Box 54806 RAK FZE , UAE

Email Us At:

sales@acculabusa.com middleeast@acculabusa.com

Phone Support:

© Copyright 2014 by ACCULAB_USA - All rights reserved